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Referenzen

- World Health Organization. Global guidelines for the prevention of surgical site infection, 2nd edition. <https://apps.who.int/iris/handle/10665/277399>. 2018.
- Suetens C, Latour K, Kärki T, et al. Prevalence of healthcare-associated infections, estimated incidence and composite antimicrobial resistance index in acute care hospitals and long-term care facilities: results from two European point prevalence surveys, 2016 to 2017. Euro Surveill 2018;23.
- World Health Organization. Decontamination and reprocessing of medical devices for health-care facilities. 2016.
- Chiu KW, Lu LS, Chiou SS. High-level disinfection of gastrointestinal endoscope reprocessing. World J Exp Med 2015;5:33-9.
- DiazGranados CA, Jones MY, Kongphet-Tran T, et al. Outbreak of Pseudomonas aeruginosa infection associated with contamination of a flexible bronchoscope. Infect Control Hosp Epidemiol 2009;30:550-5.
- Lee YK, Park JB. Steps of reprocessing and equipments. Clin Endosc 2013;46:274-9.
- McDonnell G, Burke P. Disinfection: is it time to reconsider Spaulding? J Hosp Infect 2011;78:163-70.
- Rutala WA, Weber DJ. New developments in reprocessing semicritical items. Am J Infect Control 2013;41:S60-6.
- Agency for Toxic Substances & Disease Registry. Medical Management Guidelines for Ethylene Oxide.
- World Health Organization. Concise International Chemical Assessment Document 54: Ethylene Oxide. 2003.
- Advanced Sterilization Products. Now is the Right Time to Go ETO Free.
- Centers for Disease Control and Prevention. Guideline for Disinfection and Sterilization in Healthcare Facilities, 2019.
- Muscarella LF. Use of ethylene-oxide gas sterilisation to terminate multidrug-resistant bacterial outbreaks linked to duodenoscopes. BMJ Open Gastroenterol 2019;6:e000282.
- Advanced Sterilization Products. ASP Ecosystem Brochure. Innovation to advance the Quality of Patient Care. AD-170035-01.
- Advanced Sterilization Products. STERRAD™ 100NX Sterilization System User's Guide. 2020.
- Advanced Sterilization Products. STERRAD NX™ with ALLClear™ Technology Data Sheet. AD-170009-01.
- Advanced Sterilization Products. STERRAD™ Superiority Factsheet. Volume 2 - Defining Chamber Capacity.
- Advanced Sterilization Products. STERRAD™ Systems with ALLClear™ Technology Product Brochure. AD-160029-01.
- Advanced Sterilization Products. STERRAD™ Superiority Factsheet. Volume 4 - The Importance of Instrument Validation.
- Advanced Sterilization Products. ASP ACCESS™ Technology Product Brochure. AD-160030-01.
- Advanced Sterilization Products. STERRAD VELOCITY™ Biological Indicator (BI)/Process Challenge Device (PCD) Technical Dossier.
- Advanced Sterilization Products. STERRAD VELOCITY™ Biological Indicator System Product Brochure. AD-170021-01.
- Rutala WA, Kanamori H, Sickbert-Bennett EE, et al. What's new in reprocessing endoscopes: Are we going to ensure "the needs of the patient come first" by shifting from disinfection to sterilization? Am J Infect Control 2019;47s:A62-a66.
- Culver DA, Gordon SM, Mehta AC. Infection control in the bronchoscopy suite: a review of outbreaks and guidelines for prevention. Am J Respir Crit Care Med 2003;167:1050-6.
- Hildebrand EA, Branaghan RJ, Neuman BL, et al. An Expert Perspective of Errors in Endoscope Reprocessing. Proceedings of the Human Factors and Ergonomics Society Annual Meeting 2011;55:748-752.
- Kovaleva J, Peters FTM, van der Mei HC, et al. Transmission of infection by flexible gastrointestinal endoscopy and bronchoscopy. Clinical microbiology reviews 2013;26:231-254.
- U.S. Food & Drug Administration. Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling, 2015.
- Hamasuna R, Takahashi S, Yamamoto S, et al. Guideline for the prevention of health care-associated infection in urological practice in Japan. Int J Urol 2011;18:495-502.
- Koo VSW, O'Neill P, Elves A. Multidrug-resistant NDM-1 Klebsiella outbreak and infection control in endoscopic urology. BJU International 2012;110:E922-E926.
- Advanced Sterilization Products. Let Your STERRAD™ System Do More for You.
- Advanced Sterilization Products. Cost Savings and Terminal Sterilization Among Key Benefits of the STERRAD™ Sterilization System. AD-53335-001.
- McCreanor V, Graves N. An economic analysis of the benefits of sterilizing medical instruments in low-temperature systems instead of steam. Am J Infect Control 2017;45:756-760.
- Advanced Sterilization Products. A STELLAR Study: Steam Reprocessing of Reusable Laryngoscopes and the Potential Extension of Laryngoscope Lifetimes Through a STERRAD™ Systems Alternative.
- Rogers WJ. 2 - Steam and dry heat sterilization of biomaterials and medical devices. In: Lerouge S, Simmons A, eds. Sterilisation of Biomaterials and Medical Devices: Woodhead Publishing, 2012:20-55.
- Timm D, Gonzales, D. Effect of sterilization on microstructure and function of microsurgical scissors. Surgical Services Management 1997;3:47-49.
- Advanced Sterilization Products. Assessment of Operating Costs Due to Energy and Water Use During Terminal Sterilization with STERRAD™ Systems Compared to a Steam Sterilizer.
- IARC Working Group on the Evaluation of Carcinogenic Risk to Humans. Ethylene Oxide. Chemical Agents and Related Occupations. IARC Monographs on the Evaluation of Carcinogenic Risks to Humans, No. 100F., 2012.
- Steenland K, Stayner L, Daddens J. Mortality analyses in a cohort of 18 235 ethylene oxide exposed workers: follow up extended from 1987 to 1998. Occupational and environmental medicine 2004;61:2-7.
- Ari S, Caca I, Sahin A, et al. Toxic anterior segment syndrome subsequent to pediatric cataract surgery. Cutan Ocul Toxicol 2012;31:53-7.
- Kanemitsu K, Imasaka T, Ishikawa S, et al. A comparative study of ethylene oxide gas, hydrogen peroxide gas plasma, and low-temperature steam formaldehyde sterilization. Infect Control Hosp Epidemiol 2005;26:486-9.
- Adler S, Scherrer M, Daschner FD. Costs of low-temperature plasma sterilization compared with other sterilization methods. Journal of Hospital Infection 1998;40:125-134.
- Advanced Sterilization Products. STERRAD™ Heritage Systems.
- Advanced Sterilization Products. Comparison Study of Environmental Hydrogen Peroxide Levels of STERRAD™ Systems and STERIS V-PRO® Low Temperature Sterilizers Reveals Striking Differences.
- Schneider PM. Evaluation of a new rapid readout biological indicator for use in 132°C and 135°C vacuum-assisted steam sterilization cycles. Am J Infect Control 2014;42:e17-21.
- Advanced Sterilization Products. STERRAD™ 100NX System Cycle Selection.
- STERIS. V-PRO® maX Sterilizer Cycle Overview.

asp.com

ASP Advanced Sterilization Products

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AD-210035-03-CT_B-MDR

CE
0123

ADVANCED STERILIZATION PRODUCTS, INC.
33 Technology Drive, Irvine CA 92618, USA

EC REP **ASP, The Netherlands BV**
BIC 1, 5657 BX Eindhoven, The Netherlands

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©ASP 2023. Alle Rechte vorbehalten.
ASP Deutschland GmbH