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
Quellen

1. Agency for Toxic Substances & Disease Registry. Medical Management Guidelines for Ethylene Oxide.
2. Agency for Toxic Substances & Disease Registry. Medical Management Guidelines for Formaldehyde.
3. Agency for Toxic Substances & Disease Registry. Medical Management Guidelines for Hydrogen Peroxide.
4. IARC Working Group on the Evaluation of Carcinogenic Risk to Humans. Ethylene Oxide. Chemical Agents and Related Occupations. IARC Monographs on the Evaluation of Carcinogenic Risks to Humans, No. 100F., 2012.
5. IARC Working Group on the Evaluation of Carcinogenic Risk to Humans. Formaldehyde. Chemical Agents and Related Occupations. IARC Monographs on the Evaluation of Carcinogenic Risks to Humans, No. 100F., 2012.
6. AFNOR. FD S98-135. Stérilisation des dispositifs médicaux. Guide pour la maîtrise des traitements appliqués aux dispositifs médicaux réutilisables, 2005.
7. Segurança do Paciente. Portaria Interministerial No 482, de 16 de Abril de 1999.
8. ZENCİROĞLU D. ETİLEN OKSİT STERİLİZASYONU, 2005.
9. Environmental Protection Agency. National Emission Standards for Hospital Ethylene Oxide Sterilizers. Vol 72; December 28, 2007.
10. Centers for Disease Control and Prevention. Guideline for Disinfection and Sterilization in Healthcare Facilities, 2019.
11. Advanced Sterilization Products. Voice of Customer. Urologists.
12. Kanemitsu K, Imasaka T, Ishikawa S, et al. A comparative study of ethylene oxide gas, hydrogenperoxide gas plasma, and low-temperature steamformaldehyde sterilization. Infection control and hospital epidemiology 2005;26:486-489.
13. Chobin NG. Cost analysis of three low-temperature sterilization systems at Saint Barnabas Medical Center. JHealthc Mater Manage 1994;12:29, 32-4.
14. King JS, Pink MM, Jobe CM. Assessment of reprocessed arthroscopic shaver blades. Arthroscopy: the journal of arthroscopic & related surgery: official publication of the Arthroscopy Association of North America and the International Arthroscopy Association 2006;22:1046-1052.
15. Timm D, Gonzales, D. Effect of sterilization on microstructure and function of microsurgical scissors. Surgical Services Management 1997;3:47-49.
16. Advanced Sterilization Products. Assessment of Operating Costs Due to Energy and Water Use During Terminal Sterilization with STERRAD™ Systems Compared to a Steam Sterilizer.
17. U.S. Department of Energy. Federal Energy Management Program. Water Efficiency Case Study-Veteran Affairs Hospital. 2007.
18. Steenland K, Stayner L, Deddens J. Mortality analyses in a cohort of 18 235 ethylene oxide exposed workers: follow up extended from 1987 to 1998. Occupational and environmental medicine 2004;61:2-7.
19. American Chemistry Council. Ethylene Oxide. Product Stewardship Manual - 3rd Edition.
20. Lab Manager. Safe Use of Autoclave Sterilizers, 2018.
21. Environmental Protection Agency. Evaluation of the Inhalation Carcinogenicity of Ethylene Oxide, 2016.
22. Adler S, Scherrer M, Daschner FD. Costs of low-temperature plasma sterilization compared with other sterilization methods. Journal of Hospital Infection 1998;40:125-134.
23. ACGIH®. Hydrogen Peroxide: TLV® Chemical Substances 7th Edition Documentation.
24. Advanced Sterilization Products. Comparison Study of Environmental Hydrogen Peroxide Levels of STERRAD™ Systems and STERIS V-PRO® Low Temperature Sterilizers Reveals Striking Differences.
25. World Health Organization. Ethylene Oxide.
26. Ari S, Caca I, Sahin A, et al. Toxic anterior segment syndrome subsequent to pediatric cataract surgery. Cutan Ocul Toxicol 2012;31:53-7.
27. Kanemitsu K, Kunishima H, Saga T, et al. Residual formaldehyde on plastic materials and medical equipment following low-temperature steam and formaldehyde sterilization. J Hosp Infect 2005;59:361-4.
28. Ikarashi Y, Tsuchiya T, Nakamura A. Cytotoxicity of medical materials sterilized with vapour-phase hydrogenperoxide. Biomaterials 1995;16:177-183.
29. Vink P. Residual formaldehyde in steam-formaldehyde sterilized materials. Biomaterials 1986;7:221-224.
30. Kovaleva J, Peters FTM, van der Mei HC, et al. Transmission of infection by flexible gastrointestinal endoscopy and bronchoscopy. Clinical microbiology reviews 2013;26:231-254.
31. U.S. Food and Drug Administration. Liquid Chemical Sterilization.
32. Swenson D, Wilder JA, Hancock CO. Steam sterilization validation for implementation of parametric release at a healthcare facility. Biomed Instrum Technol 2010;44:166-74.
33. McIsaac DI, Abdulla K, Yang H, et al. Association of delay of urgent or emergency surgery with mortality and use of health care resources: a propensity score-matched observational cohort study. Cmaj 2017;189:E905-e912.
34. Environmental Protection Agency. Hazardous Air Pollutants: Ethylene Oxide, 2020.
35. Joint Service Pollution Prevention Opportunity Handbook. Low Temperature Oxidative Sterilization Methods for Sterilizing Medical Devices. 2004.
36. Lesser MP. OXIDATIVE STRESS IN MARINE ENVIRONMENTS: Biochemistry and Physiological Ecology. Annual Review of Physiology 2006;68:253-278.
37. Toxics Use Reduction Institute. Formaldehyde Fact Sheet.
38. Environmental Protection Agency. Effects of Acid Rain.
39. Advanced Sterilization Products. Now is the Right Time to Go EtO Free.
40. Lerouge S, Guignot C, Tabrizian M, et al. Plasma-based sterilization: Effect on surface and bulk properties and hydrolytic stability of reprocessed polyurethane electrophysiology catheters. Journal of Biomedical Materials Research 2000;52:774-782.
41. McCreanor V, Graves N. An economic analysis of the benefits of sterilizing medical instruments in low-temperature systems instead of steam. Am J Infect Control 2017;45:756-760.
42. Advanced Sterilization Products. A STELLAR Study: Steam Reprocessing of Reusable Laryngoscopes and the Potential Extension of Laryngoscope Lifetimes Through a STERRAD™ Systems Alternative.
43. Advanced Sterilization Products. Let Your STERRAD™ System Do More for You.

AD-200077-03-CT_A-MDR 

0123

asp.com