

VOLUME **1**

STERRAD™

SUPERIORITY A FACTSHEET BY ASP™



The Power of Plasma

How do STERRAD™ with ALLClear™ Technology Sterilization Systems Operate?


- ✓ STERRAD™ with ALLClear™ Technology Sterilization Systems utilise a combination of **hydrogen peroxide (H₂O₂)** and **low-temperature gas plasma** to rapidly and safely sterilize validated medical devices and materials, without leaving toxic residues.¹
- ✓ This approach offers **safety and efficiency advantages** over alternative reprocessing modalities such as ethylene oxide (EtO), formaldehyde gas (FO) and other H₂O₂ systems, as described below.



ASP pioneered the use of gas plasma in low temperature sterilizers in the 1990s and has led the market since then.”

How Does Plasma Work?

1. H₂O₂, combined with low-temperature **gas plasma**, results in the production of microbicidal free radicals, such as hydroxyl, which **disrupt essential cell components**
2. Plasma dissociates unreacted H₂O₂ into oxygen and water and **removes any residual H₂O₂** from the load



Safe for the User

Commonly used sterilants, including EtO, FO and H₂O₂, and their residues are associated with acute and long-term toxicities, such as chronic irritation, central nervous system depression, severe allergic reactions and many more.^{2,3,4,5,6}

→ Sterilizers which lack gas plasma technology **risk exposing users to such residues** or require **lengthy and disruptive measures** to eliminate them


→ Sterilization with **EtO requires aeration** of the load to help keep staff safe, with typical sterilization cycles lasting **16–17 hours**⁷

→ Sterilization with H₂O₂ does not typically require aeration of the load before handling, however, sterilizer models without gas plasma technology have been shown to produce **H₂O₂ emissions above that deemed safe** by the American Conference of Governmental Industrial Hygienists (ACGIH®)⁸


In contrast, by utilising gas plasma technology to remove residual H₂O₂, STERRAD™ with ALLClear™ Technology Sterilization Systems **reduce exposure to harmful residues, to safe levels**

In line with this, STERRAD™ with ALLClear™ Technology Sterilization Systems, H₂O₂ emissions, at the user's breathing zone level, are up to 67 times lower than for STERIS V-PRO® sterilizers:⁹


1. STERRAD™ 100NX with ALLClear™ Technology System FLEX Cycles




2. V-PRO® maX Flexible Cycles



3. STERRAD™ 100NX with ALLClear™ Technology System STANDARD Cycles



4. V-PRO® maX lumen Cycles



Key: ■ Above ACGIH® Permissible limits (5ppm) ■ Peak concentration (ppm)



STERIS V-PRO® sterilizers showed instantaneous peak measurements of H₂O₂ reached up to 20ppm at the user's breathing zone level, contributing to a more hazardous working environment⁹



Safe for the Patient

- ✓ Sterilizing agents, including EtO, FO and H₂O₂, can leave toxic residues on medical devices,^{10,14} putting patients at risk of health complications^{10,11}
- ✓ By utilising gas plasma technology, STERRAD™ with ALLClear™ Technology Sterilization Systems leave medical devices free from toxic residues¹, ensuring that they are safe for the patient
- ✓ In addition, lengthy instrument turnover times, such as those associated with EtO, may delay operating schedules due to unavailability of surgical instruments, compromising patient safety¹⁵



1-2% EtO¹⁰

Concentrations of unchanged EtO have been measured in sterilized devices.



19/893 eyes had TASS¹¹

TASS has occurred as a result of EtO-sterilized vitrectomy packs. No TASS cases were observed with non-EtO sterilized packs.

EtO and H₂O₂ can have a toxic impact on the environment and so need to be strictly regulated in order to reduce their environmental impact^{16,17}
 Gas plasma eliminates potentially harmful emissions, allowing efficiency without compromising the environment

Safe for the Environment



Key TAKEAWAYS

- ✓ STERRAD™ Sterilization systems with ALLClear™ Technology minimize exposure to harmful sterilizing agents residues by utilizing plasma.
- ✓ Other sterilizers, which utilize H₂O₂ without gas plasma technology, produce H₂O₂ emissions that exceed recommended exposure thresholds. STERRAD™ Sterilization systems with ALLClear™ Technology reduce these emissions to safe levels, ensuring that users, patients and the environment are protected, without the need of weekly chamber leak tests or quarterly gasket inspection required.

asp.com

ASP Advanced Sterilization Products

ASP International GmbH, Zug Branch
 Bahnhofstrasse 2, Zug 6300, Switzerland
 ©ASP 2022. All Rights Reserved.



ADVANCED STERILIZATION PRODUCTS, INC.
 33 Technology Drive, Irvine CA 92618, USA



ASP, The Netherlands BV
 BIC 1, 5657 BX Eindhoven, The Netherlands

CE
 0123

1. Advanced Sterilization Products. STERRAD 100NX™ User Guide (page 8). 2. Agency for Toxic Substances & Disease Registry, Medical Management Guidelines for Ethylene Oxide. 3. IARC Working Group on the Evaluation of Carcinogenic Risk to Humans. Ethylene Oxide. Chemical Agents and Related Occupations. IARC Monographs on the Evaluation of Carcinogenic Risks to Humans, No. 100F, 2012. 4. Agency for Toxic Substances & Disease Registry, Medical Management Guidelines for Formaldehyde. 5. IARC Working Group on the Evaluation of Carcinogenic Risk to Humans. Formaldehyde, Chemical Agents and Related Occupations. IARC Monographs on the Evaluation of Carcinogenic Risks to Humans, No. 100F, 2012. 6. Agency for Toxic Substances & Disease Registry, Medical Management Guidelines for Hydrogen Peroxide. 7. Kanemitsu K, Imasaka T, Ishikawa S, et al. A comparative study of ethylene oxide gas, hydrogen peroxide gas plasma, and low-temperature steam formaldehyde sterilization. Infection control and hospital epidemiology 2005;26:486-489. 8. ACGIH®. Hydrogen Peroxide: TLY™ Chemical Substances 7th Edition Documentation. 9. Advanced Sterilization Products. Comparison Study of Environmental Hydrogen Peroxide Levels of STERRAD™ Systems and STERIS V-PRO® Low Temperature Sterilizers Reveals Striking Differences. 10. World Health Organization. Ethylene Oxide. 11. Ari S, Caca I, Sahin A, et al. Toxic anterior segment syndrome subsequent to pediatric cataract surgery. Cutan Ocul Toxicol 2012;31:53-7. 12. Kanemitsu K, Kunishima H, Saga T, et al. Residual formaldehyde on plastic materials and medical equipment following low-temperature steam and formaldehyde sterilization. J Hosp Infect 2005;59:361-4. 13. Vink P. Residual formaldehyde in steam-formaldehyde sterilized materials. Biomaterials 1986;7:221-224. 14. Ikarashi Y, Tsuchiya T, Nakamura A. Cytotoxicity of medical materials sterilized with vapour-phase hydrogen peroxide. Biomaterials 1995;16:177-183. 15. Mclsaac DI, Abdulla K, Yang H, et al. Association of delay of urgent or emergency surgery with mortality and use of health care resources: a propensity score-matched observational cohort study. Cmaj 2017;189:E905-e912. 16. Lesser MR. OXIDATIVE STRESS IN MARINE ENVIRONMENTS: Biochemistry and Physiological Ecology. Annual Review of Physiology 2006;68:253-278. 17. Environmental Protection Agency. Hazardous Air Pollutants: Ethylene Oxide, 2020.