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Références

- Chobin NG. Cost analysis of three low-temperature sterilization systems at Saint Barnabas Medical Center. *J Healthc Mater Manage.* 1994;12:29, 32-4.
- Swenson D, Wilder JA, Hancock CO. Steam sterilization validation for implementation of parametric release at a healthcare facility. *Biomed Instrum Technol.* 2010;44:166-74.
- Kanemitsu K, Imasaka T, Ishikawa S, et al. A comparative study of ethylene oxide gas, hydrogen peroxide gas plasma, and low-temperature steam formaldehyde sterilization. *Infect Control Hosp Epidemiol.* 2005;26:486-9.
- McIsaac DI, Abdulla K, Yang H, et al. Association of delay of urgent or emergency surgery with mortality and use of health care resources: a propensity score-matched observational cohort study. *Cmaj.* 2017;189:E905-e912.
- Teillant A, Gandra S, Barter D, et al. Potential burden of antibiotic resistance on surgery and cancer chemotherapy antibiotic prophylaxis in the USA: a literature review and modelling study. *Lancet Infect Dis.* 2015;15:1429-37.
- Agency for Healthcare Research and Quality. 2013 Annual Hospital-Acquired Condition Rate and Estimates of Cost Savings and Deaths Averted From 2010 to 2013. Last updated [Accessed online at https://www.ahrq.gov/sites/default/files/publications/files/hacrate2013_0.pdf on 13/10/2020].
- Dancer SJ, Stewart M, Coulombe C, et al. Surgical site infections linked to contaminated surgical instruments. *J Hosp Infect.* 2012;81:231-8.
- Tosh PK, Disbot M, Duffy JM, et al. Outbreak of *Pseudomonas aeruginosa* surgical site infections after arthroscopic procedures: Texas, 2009. *Infect Control Hosp Epidemiol.* 2011;32:1179-86.
- Kirkland KB, Briggs JP, Trivette SL, et al. The impact of surgical-site infections in the 1990s: attributable mortality, excess length of hospitalization, and extra costs. *Infect Control Hosp Epidemiol.* 1999;20:725-30.
- Guédon ACP, Wauben LSG, van der Eijk AC, et al. Where are my instruments? Hazards in delivery of surgical instruments. *Surgical endoscopy.* 2016;30:2728-2735.
- Timm D, Gonzales, D. Effect of sterilization on microstructure and function of microsurgical scissors. *Surgical Services Management.* 1997;3:47-49.
- Kovaleva J, Peters FTM, van der Mei HC, et al. Transmission of infection by flexible gastrointestinal endoscopy and bronchoscopy. *Clinical microbiology reviews.* 2013;26:231-254.
- Advanced Sterilization Products. Assessment of Operating Costs Due to Energy and Water Use During Terminal Sterilization with STERRAD™ Systems Compared to a Steam Sterilizer.
- U.S. Department of Energy. Federal Energy Management Program. Water Efficiency Case Study - Veteran Affairs Hospital. 2007.
- Agency for Toxic Substances & Disease Registry. Medical Management Guidelines for Ethylene Oxide. Last updated [Accessed online at <https://www.atsdr.cdc.gov/mmg/mmg.asp?id=730&tid=133> on 30/06/2020].
- Agency for Toxic Substances & Disease Registry. Medical Management Guidelines for Formaldehyde. Last updated [Accessed online at <https://www.atsdr.cdc.gov/mmg/mmg.asp?id=216&tid=39> on 30/06/2020].
- Agency for Toxic Substances & Disease Registry. Medical Management Guidelines for Hydrogen Peroxide. Last updated [Accessed online at <https://www.atsdr.cdc.gov/mmg/mmg.asp?id=304&tid=55> on 30/06/2020].
- IARC Working Group on the Evaluation of Carcinogenic Risk to Humans. Ethylene Oxide. Chemical Agents and Related Occupations. IARC Monographs on the Evaluation of Carcinogenic Risks to Humans, No. 100F., 2012.
- IARC Working Group on the Evaluation of Carcinogenic Risk to Humans. Formaldehyde. Chemical Agents and Related Occupations. IARC Monographs on the Evaluation of Carcinogenic Risks to Humans, No. 100F., 2012.
- AFNOR. FD S98-135. Stérilisation des dispositifs médicaux. Guide pour la maîtrise des traitements appliqués aux dispositifs médicaux réutilisables. Last updated [Accessed online at <https://www.boutique.afnor.org/norme/fd-s98-135/sterilisation-des-dispositifs-medicaux-guide-pour-la-maitrise-des-traitements-appliques-aux-dispositifs-medicaux-reutilisables/article/776199/fa136818> on 30/07/2020].
- Segurança do Paciente. Portaria Interministerial No 482, de 16 de Abril de 1999. Last updated [Accessed online at <https://www20.anvisa.gov.br/segurancadopaciente/index.php/legislacao/item/portaria-interministerial-n-482-de-16-de-abril-de-1999> on 30/07/2020].
- ZENCİROĞLU D. ETİLEN OKSİT STERİLİZASYONU. Last updated [Accessed online at <https://www.das.org.tr/kitaplar/kitap2005/11-05.pdf> on 30/07/2020].
- Environmental Protection Agency. National Emission Standards for Hospital Ethylene Oxide Sterilizers. Vol 72; December 28, 2007. Last updated [Accessed online at <https://www3.epa.gov/ttn/atw/area/fr28de07b.pdf> on 30/07/2020].
- Centers for Disease Control and Prevention. Guideline for Disinfection and Sterilization in Healthcare Facilities. Last updated [Accessed online at <https://www.cdc.gov/infectioncontrol/pdf/guidelines/disinfection-guidelines-H.pdf> on 30/07/2020].
- Adler S, Scherrer M, Daschner FD. Costs of low-temperature plasma sterilization compared with other sterilization methods. *Journal of Hospital Infection.* 1998;40:125-134.
- Advanced Sterilization Products. Comparison Study of Environmental Hydrogen Peroxide Levels of STERRAD™ Systems and STERIS V-PRO™ Low Temperature Sterilizers Reveals Striking Differences. Last updated [Accessed online at [https://www.asp.com/sites/default/files/pdf/054380-180619_STERRAD_H2O2_Case_Study%20\(3\).pdf](https://www.asp.com/sites/default/files/pdf/054380-180619_STERRAD_H2O2_Case_Study%20(3).pdf) on 30/06/2020].
- World Health Organization. Report on the Burden of Endemic Health Care-Associated Infection Worldwide. Last updated [Accessed online at https://apps.who.int/iris/bitstream/handle/10665/80135/9789241501507_eng.pdf;jsessionid=86B33E75BBE594E6C2CBD998F1C88CA5?sequence=1, on 13/10/2020].
- World Health Organization. Health care-associated infections FACT SHEET. Last updated [Accessed online at https://www.who.int/gpsc/country_work/gpsc_ccisc_fact_sheet_en.pdf?ua=1 on 13/10/2020].
- Emerson CB, Eyzaguirre LM, Albrecht JS, et al. Healthcare-associated infection and hospital readmission. *Infection control and hospital epidemiology.* 2012;33:539-544.
- Lamarsalle L, Hunt B, Schauf M, et al. Evaluating the clinical and economic burden of healthcare-associated infections during hospitalization for surgery in France. *Epidemiol Infect.* 2013;141:2473-82.
- World Health Organization. Guidelines on Core Components of Infection Prevention and Control Programmes at the National and Acute Health Care Facility Level. Last updated [Accessed online at https://www.ncbi.nlm.nih.gov/books/NBK401773/pdf/Bookshelf_NBK401773.pdf on 13/10/2020].
- European Centre for Disease Prevention and Control. Annual Epidemiological Report 2016 - Healthcare-associated infections acquired in intensive care units. Last updated [Accessed online at https://www.ecdc.europa.eu/sites/portal/files/documents/AER-HCAI_ICU_3_0.pdf on 13/10/2020].
- Schneider PM. Evaluation of a new rapid readout biological indicator for use in 132°C and 135°C vacuum-assisted steam sterilization cycles. *Am J Infect Control.* 2014;42:e17-21.
- Advanced Sterilization Products. Cost Savings and Terminal Sterilization Among Key Benefits of the STERRAD™ Sterilization System. AD-53335-001.
- McCreaon V and Graves N. An economic analysis of the benefits of sterilizing medical instruments in low-temperature systems instead of steam. *Am J Infect Control.* 2017;45:756-760.
- Advanced Sterilization Products. A STELLAR Study: Steam Reprocessing of Reusable Laryngoscopes and the Potential Extension of Laryngoscope Lifetimes Through a STERRAD™ Systems Alternative.
- Advanced Sterilization Products. Let Your STERRAD™ System Do More for You.
- Graves N, Weinhold D, Tong E, et al. Effect of healthcare-acquired infection on length of hospital stay and cost. *Infect Control Hosp Epidemiol.* 2007;28:280-92.
- Nursing in Practice. Healthcare-associated infections in the community. Last updated [Accessed online at <https://www.nursinginpractice.com/article/healthcare-associated-infections-community> on 13/10/2020].
- CareFusion. Taking a zero tolerance approach to preventable surgical site infections in UK hospitals. Last updated [Accessed online at http://www.carefusion.co.uk/documents/international/continuing-education/infection-prevention/IP_Under-the-Knife_CE_EN.pdf on 13/10/2020].

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ADVANCED STERILIZATION PRODUCTS, INC.
33 Technology Drive, Irvine CA 92618, USA

ASP, The Netherlands BV
BIC 1, 5657 BX Eindhoven, The Netherlands

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©ASP 2022. Tous droits réservés.
ASP France SAS