CIDEX® OPA Solution Efficacy Against Human Coronavirus

Background

Coronaviruses are a family of enveloped positive-sense single-stranded RNA viruses, capable of rapid mutation and recombination. Coronaviruses are a common cause of mild upper respiratory infections, but they also include more serious illness like SARS and MERS, which can lead to severe disease. Coronaviruses are common in many mammals and, although rare, animal coronaviruses can infect humans and then spread between people¹.

The 2019 novel coronavirus (2019-nCoV) is a new virus that has confirmed human-to-human transmission². It was first identified during an outbreak in Wuhan, Hubei Province, China in late 2019³. Although most patients with coronavirus have mild symptoms, the 2019-nCoV has resulted in severe cases sometimes escalating to pneumonia. Deaths have also been reported within China⁴.

The exact way 2019-nCoV is spread is not yet fully understood. It is suspected that it is spread in a similar manner as other virulent strains of coronaviruses, such as MERS and SARS, via respiratory droplets produced when an infected person coughs or sneezes. There also may be some spread when a person touches a surface or object that has virus on it and then touches his or her own mouth, nose, or eyes. Spread of SARS and MERS between people has generally occurred between close contacts⁴. SARS and MERS viruses can also survive on hard surfaces for greater than 48 hours⁵. It is likely that 2019-nCoV can survive in a similar manner on hard surfaces.

Efficacy against Coronavirus

Although ASP has not specifically conducted testing of CIDEX OPA Solution against 2019-nCoV, CIDEX OPA Solution has been tested for efficacy against human coronavirus. Product development testing was performed using the United States Environmental Protection Agency (EPA) Virucide Assay Method⁶. In this development testing, preparations of viral samples of human coronavirus were inactivated by CIDEX OPA Solution after a 5-minute exposure time at 20°C; and exceeded the required $\geq 3\text{-log}_{10}$ reduction in virus titer. Testing was also completed using a diluted strength CIDEX OPA formulation to present adequate challenge to CIDEX OPA. The results show that even at a diluted concentration, CIDEX OPA Solution inactivates Human Coronavirus. These test results satisfy EPA requirements.

Based on this testing data regarding the efficacy of CIDEX OPA against coronaviruses, 2019-nCoV should be susceptible to CIDEX OPA.

Additionally, coronaviruses such as 2019-nCoV are lipid viruses (enveloped viruses). Lipid viruses are typically susceptible to various low and medium disinfection modalities according to the guide of the hierarchy of biocide resistance published by the Centers for Disease Control and Prevention (CDC)⁷. CIDEX OPA is a High-Level-Disinfectant and has been shown to be effective against microorganisms such as mycobacteria and polio virus, which are microorganisms that typically exhibit greater resistance to disinfection modalities than lipid viruses like the coronavirus. Based on what is currently known about microorganism resistance to disinfection, CIDEX OPA would be efficacious against lipid viruses such as coronavirus.

Finally, it is important to properly follow the CIDEX OPA solution Instruction for Use (IFU) for detailed directions when processing any instrument using CIDEX OPA, especially the pre-cleaning of all instruments as well as rinsing procedures within the IFU to ensure complete high-level disinfection. HLD guidelines in your facilities as well as the HLD instructions of manufacturers should also be followed.

Please keep in mind to follow the safety precautions of the local health authorities of your country as well as health care facility when dealing with infected patients.

For any additional questions regarding 2019-nCoV or Coronaviruses in relation to our products, please submit a Medical Information Request by visiting www.asp.com.

Sincerely,

ASP Medical Affairs
Advanced Sterilization Products

www.ASP.com

33 Technology Drive.

Irvine, CA 92618

Email: Medical@ASP.com

Please see attached Package Insert for more information about CIDEX OPA

REFERENCES:

1. Coronavirus. <https://www.who.int/health-topics/coronavirus>. Accessed 01/28/20.
2. Chan JF-W, Yuan S, Kok K-H, et al. *A familial cluster of pneumonia associated with the 2019 novel coronavirus indicating person-to-person transmission: a study of a family cluster*. Lancet. 2020; (published online Jan 24). [https://doi.org/10.1016/S0140-6736\(20\)30154-9](https://doi.org/10.1016/S0140-6736(20)30154-9)
3. About 2019 Novel Coronavirus (2019-nCoV). <https://www.cdc.gov/coronavirus/2019-ncov/about/index.html>. Accessed 01/28/2020.
4. Chen Wang, Peter W Horby, Frederick G Hayden, George F Gao. *A novel coronavirus outbreak of global health concern*. Lancet. 2020; (published online Jan 24). [https://doi.org/10.1016/S0140-6736\(20\)30185-9](https://doi.org/10.1016/S0140-6736(20)30185-9).
5. Otter JA, Donskey C, Yezli S, Douthwaite S, Goldenberg SD, Weber DJ. Transmission of SARS and MERS coronaviruses and influenza virus in healthcare settings: the possible role of dry surface contamination. J Hosp Infect. 2016;92:235–250. doi: 10.1016/j.jhin.2015.08.027.
6. EPA Pesticide Assessment Guidelines: Subdivision G, Product Performance.
7. Decreasing order of resistance of microorganisms to disinfection and sterilization and the level of disinfection or sterilization. <https://www.cdc.gov/infectioncontrol/guidelines/disinfection/tables/figure1.html>. Accessed 01/28/20.

DISCLAIMER

Please note that medical devices or medicinal products are only approved for, and must be used in accordance with, those indications and dosages that appear in their Instruction for Use, Summary of Product Characteristics and/or Package Insert. Some treatment modalities included in the literature are the opinions of the authors only and not recommendations or guidance for use.

ADVERSE EVENT REPORTING

Remember to contact our local representative in case any adverse event or product quality complaint occur while using our product and comply with your local regulations in terms of notification.

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INSTRUCTIONS FOR USE

CIDEX® OPA *ortho*-Phthalaldehyde Solution High Level Disinfectant

Active Ingredient <i>ortho</i> -phthalaldehyde	0.55%
Inert Ingredients	99.45%
Dipotassium hydrogen phosphate	
Potassium dihydrogen phosphate	
Benzotriazole	
Citric acid	
D&C Green Dye #5	
N-(hydroxyethyl)-ethylenediaminetriacetic acid (HEDTA)	
Water	
Total	100.00%

Does not require activation before use.

INDICATIONS FOR USE

CIDEX® OPA Solution is a high level disinfectant for reprocessing heat sensitive medical devices when used according to the Directions for Use. CIDEX® OPA Solution is intended for use in manual (bucket and tray) systems made from polypropylene, acrylonitrile-butadiene-styrene (ABS), polyethylene, glass-filled polypropylene and/or polycarbonate plastics. CIDEX® OPA Solution may also be used in automated endoscope reprocessors according to the manufacturer's instructions and should be monitored with CIDEX® OPA Solution Test Strips. See DIRECTIONS FOR USE – Reusage for Disinfection.

Medical devices reprocessed in CIDEX® OPA Solution must first be cleaned according to a validated cleaning protocol or standard.

Antimicrobial Level of Activity: CIDEX® OPA Solution can be used at the following antimicrobial level of activity:

High Level Disinfectant: CIDEX® OPA Solution is a high level disinfectant for semi-critical devices when used or reused, according to the Directions for Use, at or above its Minimum Effective Concentration (MEC) of 0.3% as determined by CIDEX® OPA Solution Test Strips, at 20°C (68°F) with an immersion time of at least 5 minutes for a reuse period of up to 14 days.

Reuse Period for Disinfection: CIDEX® OPA Solution has demonstrated disinfection efficacy in the presence of 5% organic soil contamination and microbiological burden during reuse. CIDEX® OPA Solution may be reused for up to a maximum of 14 days provided

the required conditions of *ortho*-phthalaldehyde concentration and temperature exist based upon monitoring described in the Directions for Use. DO NOT rely solely on days in use. Concentration of this product during its reuse life must be verified by the CIDEX® OPA Solution Test Strip prior to each use to determine that the concentration of *ortho*-phthalaldehyde is above the MEC of 0.3%. The product must be discarded after 14 days, even if the CIDEX® OPA Solution Test Strip indicates a concentration above the MEC.

General Information on Selection and Use of Disinfectants for Medical Device Reprocessing: Choose a disinfectant with the level of antimicrobial activity that is appropriate for the reusable device. Follow the reusable device labeling and standard institutional practices. In the absence of complete instructions, use the following process:

First, for patient contacting devices, determine whether the reusable device to be reprocessed is a critical or semi-critical device.

Critical device: *Presents a high risk of infection if not sterile. Routinely penetrates the skin or mucous membranes during use or is otherwise used in normally sterile tissue of the body.*

Semi-critical device: *Makes contact with mucous membranes but does not ordinarily penetrate normally sterile areas of the body.*

Second, determine if sterilization or high level disinfection is required.

Critical device (e.g. laparoscopes and microsurgical instruments): *Sterilization is required.*

Semi-critical reusable device (e.g. gastroscopes): *Sterilization is required whenever feasible; where not feasible, high level disinfection is the minimum acceptable process.*

Third, select a disinfectant that is labeled for the appropriate antimicrobial level and is compatible with the reusable device. Follow directions for the disinfectant.

Microbiocidal Activity: The following table indicates the spectrum of activity as demonstrated by testing of CIDEX® OPA Solution using prescribed test methods.

MICROORGANISM

VEGETATIVE MICROORGANISMS

Staphylococcus aureus

Salmonella enterica (formerly Salmonella choleraesuis)

Pseudomonas aeruginosa

Mycobacterium bovis

FUNGI

Trichophyton mentagrophytes

VIRUSES NON-ENVELOPED

Poliovirus Type 1

Rhinovirus Type 42

Adenovirus Type 2

Vaccinia (Wyeth)

Coxsackievirus Type B-3

VIRUSES ENVELOPED

Coronavirus

Cytomegalovirus

Influenza Virus [Hong Kong]

HIV-1

Herpes simplex Types 1,2

Material Compatibility: CIDEX® OPA Solution has been tested and found to be compatible with the materials shown below.

METALS¹

Aluminum

Anodized aluminum²

Brass

Carbon steel

Chrome plated brass²

Chrome plated steel²

Copper

Nickel plated brass²

Nickel silver alloy²

Stainless steel³

Titanium

Tungsten carbide²

Vanadium steel⁴

PLASTICS⁵

Polymethylmethacrylate (Acrylic)

Nylon

Polyethylene terephthalate (Polyester)

Polystyrene

Polyvinylchloride (PVC)⁶

Acrylonitrile/butadiene/styrene (ABS)

Polysulfone

Polycarbonate⁷

Polyethylene

Polypropylene

Acetal

PTFE

Polyamide

ELASTOMERS⁵

Polychloroprene (Neoprene)

Kraton G

Polyurethane

Silicone rubber⁶

Natural rubber latex

ADHESIVES⁵

Cyanoacrylate⁸

EPO-TEK 301 epoxy⁸

EPO-TEK 353 epoxy

DENTAL MATERIALS⁹

Polysulphide

Addition silicone

Polyether

1. Exposed to 31 days (744 hours) of continuous contact with CIDEX® OPA Solution with no effect unless otherwise noted.
2. Shows signs of surface discoloration at 7 days or greater.
3. Most grades tested show no effect. Others may exhibit slight discoloration at 7 days or greater. Stainless steel 440 shows rust at 14 days immersion.
4. Treated with 500 cycles of CIDEX® OPA Solution. Surface breakdown noted after 150 cycles (25 hour total contact).
5. Exposed to 7 days of continuous contact with CIDEX® OPA Solution with no effect unless otherwise noted.
6. Some grades or applications exhibit discoloration.
7. Some sonic welded parts may exhibit crazing.
8. Some loss in shear strength but show no signs of severe degradation.
9. Dental materials treated with 1 cycle of CIDEX® OPA Solution for 10 minutes with no effects.

Fujinon, Olympus and Pentax endoscopes are compatible with CIDEX® OPA Solution. If questions arise regarding the compatibility of a device with CIDEX® OPA Solution, contact the device manufacturer.

Cleaning Agent Compatibility: CIDEX® OPA Solution is compatible with enzymatic detergents which are near neutral (6-8) in pH, low foaming, and easily rinsed from equipment (e.g., CIDEZYME® and CIDEZYME® LF Enzymatic Detergents). Detergents that are either highly acidic or alkaline are not recommended as cleaning agents.

CONTRAINDICATIONS

1. CIDEX® OPA Solution should not be utilized to process any urological instrumentation used to treat patients with a history of bladder cancer. In rare instances CIDEX® OPA Solution has been associated with anaphylaxis-like reactions in bladder cancer patients under going repeated cystoscopies.
2. CIDEX® OPA Solution should not be utilized to process instrumentation for patients with known sensitivity to CIDEX® OPA Solution or any of its components.
3. CIDEX® OPA Solution should not be used to sterilize heat sensitive medical devices.

WARNINGS

CAUTION:

Contains 0.55% *ortho*-Phthalaldehyde.

Contains *ortho*-phthalaldehyde. May produce an allergic reaction, including anaphylaxis.

Safety data sheet available on request.

Read label before use.

1. May elicit an allergic reaction. Possible allergic reactions, including anaphylaxis have been reported in rare instances. In the majority of these instances health care workers were not using the product in a well-ventilated room or not wearing proper personal protective equipment. [See PRECAUTIONS].
2. Avoid contact with eyes, skin, or clothing. [See PRECAUTIONS – for important information on how to protect eyes, skin and clothing.] Direct contact with eyes may cause irritation. Direct contact with skin may cause temporary staining. Repeated contact may cause skin and respiratory sensitization. Sensitization can result in anaphylaxis in rare instances. In case of eye contact, immediately flush eyes with large quantities of water for at least 15 minutes. Seek medical attention. In case of skin contact, immediately wash with soap and water. Refer to the MSDS for additional information. Do not form sprays, mists or aerosols of this product.
3. Ingestion may cause irritation or chemical burns of the mouth, throat, esophagus and stomach. If swallowed, DO NOT INDUCE VOMITING. Drink large quantities of water and call a physician immediately. Probable mucosal damage from oral exposure may contraindicate the use of gastric lavage.
4. Avoid exposure to *ortho*-Phthalaldehyde vapors, as they may be irritating to the respiratory tract and eyes. May cause stinging sensation in the nose and throat, discharge, coughing, chest discomfort and tightness, difficulty with breathing, wheezing, tightening of throat, urticaria (hives), and/or sensitization/anaphylaxis in rare instances, rash, loss of smell, tingling of mouth or lips, dry mouth or headache. May aggravate a pre-existing asthma or bronchitis condition. In case of adverse reactions from inhalation of vapor, move to fresh air. If breathing is difficult, oxygen may be given by qualified personnel. If symptoms persist, seek medical attention.
5. The use of CIDEX® OPA Solution with semi-critical devices must be part of a validated rinsing procedure as provided by the device manufacturer. See DIRECTIONS FOR USE Rinsing Instructions for important information on rinsing. Improper rinsing can result in anaphylaxis in rare instances.
6. ALWAYS follow the Directions For Use Rinsing Instructions (Part B) and the SPECIAL INSTRUCTIONS for transesophageal echocardiography (TEE) probes in Part C EXACTLY or residues of CIDEX® OPA may remain on the device. Failure to follow rinsing instructions exactly has resulted in reports of chemical burns, irritation, and staining of the mouth, throat, oesophagus and stomach. Improper rinsing can result in anaphylaxis in rare instances.

PRECAUTIONS

Follow Hospital policy and protocol when handling and cleaning soiled devices.

1. When disinfecting devices, use gloves of appropriate type and length. Nitrile rubber gloves or butyl rubber gloves may be used. When using latex rubber gloves, the user should double glove and/or change single gloves frequently, e.g., after 10 minutes of exposure. For those individuals who are sensitive to latex or other components in latex gloves, 100% synthetic copolymer gloves are recommended.
Note: Contact with CIDEX® OPA Solution may stain exposed skin or clothing.
2. Eye protection, such as goggles or face shields, should always be used when handling CIDEX® OPA.
3. Fluid repellent gown or apron should always be used to prevent OPA from getting on clothing or skin.
4. Use CIDEX® OPA Solution in a dedicated room with adequate ventilation and in closed containers with tight-fitting lids. Please reference local/state/country requirements for appropriate number of air exchanges per hour. If adequate ventilation is not provided

- by the existing air handling system, use in local exhaust hoods, or in ductless fume hoods/portable ventilation devices which contain filter media which absorb *ortho*-phthalaldehyde from the air
- Contaminated reusable devices **MUST BE THOROUGHLY CLEANED** prior to disinfection, since residual contamination with soil or lubricants will decrease the effectiveness of the disinfectant.
 - The user **MUST** adhere to the Directions for Use, as modification to the Directions for Use may affect the safety and effectiveness of the disinfectant.
 - Do not use CIDEX® OPA Solution on medical devices that are intended for use in a sterile area of the body (e.g. cataract surgical instruments).
 - The reusable device manufacturer should provide the user with a validated reprocessing procedure for that device using CIDEX® OPA Solution.
 - The use of CIDEX® OPA Solution in automated endoscope reprocessors must be part of a validated reprocessing procedure. Use CIDEX® OPA Solution Test Strips to detect *ortho*-Phthalaldehyde concentration before each cycle to detect the MEC. Follow the Directions for Use provided with the CIDEX® OPA Solution Test Strips.
 - If using filtering face pieces, please ensure that they are capable of preventing inhalation of vapor. Do not use filtering face pieces that are only designed to prevent inhalation of particulates (e.g. do not use N95 filtering face piece respirator).
-

DIRECTIONS FOR USE

Cleaning/Decontamination: Blood, other body fluids, and lubricants must be thoroughly cleaned from the surfaces and lumens of medical devices before reprocessing in the disinfectant. Blood and other body fluids should be disposed of according to all applicable regulations for infectious waste disposal.

Refer to the reusable device manufacturer's labeling for instructions on disassembly, decontamination, cleaning and leak testing of their equipment.

Before immersion in CIDEX® OPA Solution, thoroughly clean devices, including all lumens, using a cleaning protocol or standard.

Thoroughly rinse and rough dry all surfaces and lumens of cleaned devices.

Usage: **NO ACTIVATION IS REQUIRED.**

Record the date the container was opened on the container label, or in a logbook. After opening, the solution remaining in the container may be stored for up to 75 days (providing the 75 days does not extend past the expiration date on the container) until used. Record the date the solution was poured out of the original container into a secondary container in a log book (separate from the one mentioned above), or on a label affixed to the secondary container. The solution in the secondary container can be used for a period up to 14 days. The product must be discarded after 14 days even if the CIDEX® OPA Solution Test Strip indicates a concentration above the MEC.

A. High Level Disinfection: Immerse device completely, filling all lumens and eliminating air pockets, in CIDEX® OPA Solution for a minimum of 5 minutes at 20°C (68°F) or higher to destroy all pathogenic microorganisms, including *Mycobacterium bovis*, *Pseudomonas aeruginosa*, pathogenic fungi, and viruses (Poliovirus Type 1; Adenovirus Type 2; Herpes simplex Types 1,2; HIV-1;

Influenza Type A [Hong Kong]; Vaccinia; Coronavirus; Coxsackievirus Type B-3; Cytomegalovirus; Rhinovirus Type 42). Remove device from the solution and rinse thoroughly following the rinsing instructions below.

B. Rinsing Instructions

1. RINSING PROCEDURE

a) Manual Processing:

- Following removal from CIDEX® OPA Solution, thoroughly rinse the medical device by immersing it completely in a large volume (e.g. 8 litres) of water. Use sterile water unless potable water is acceptable. See item 2 or 3 below.
- Keep the device totally immersed for a minimum of 1 minute in duration, unless a longer time is specified by the reusable device manufacturer.
- Manually flush all lumens with large volumes (not less than 100ml) of rinse water unless otherwise noted by the device manufacturer.
- Remove the device and discard the rinse water. Always use fresh volumes of water for each rinse. Do not reuse the water for rinsing or any other purpose.
- Repeat the procedure TWO (2) additional times, for a total of THREE (3) RINSES, with large volumes of fresh water to remove CIDEX® OPA Solution residues. Residues may cause serious side effects, including anaphylaxis in rare instances. SEE WARNINGS. THREE (3) SEPARATE, LARGE VOLUME WATER IMMERSION RINSES ARE REQUIRED.
- Refer to the reusable medical device manufacturer's labelling for additional rinsing instructions.

b) Automated Processing:

- Select a rinse cycle on an automatic endoscope reprocessor that has been validated for use with this product.
- Ensure that the automated rinse cycle selected will thoroughly rinse the medical device including all lumens with large volumes of sterile or potable water equivalent to the reusable device manufacturer's recommendations.
- Verify that each rinse is a minimum of 1 minute in duration unless the reusable device manufacturer specifies a longer time. Ensure that a fresh volume of water is used for each rinse. Do not reuse the water for rinsing or any other purpose.
- Refer to the reusable device manufacturer's labelling for additional rinsing instructions.

2. STERILE WATER RINSE: There are certain devices and procedures that may require sterile water rinsing after high level disinfection. The following devices should be rinsed with sterile water using clean technique:

- Devices intended for use in normally sterile areas of the body.
- Devices intended for use in known immune compromised patients, or potentially immune compromised patients based on institutional procedures (e.g., high risk population served).
- Bronchoscopes, and other semi critical devices may be at risk of contamination when rinsing with potable water. These devices should undergo a sterile water rinse. There may be microorganisms present in the potable water that may be pathogenic for immune compromised patients which may place them at high risk for infection.
- Other semi-critical, lumened instruments should be rinsed with sterile water using a sterile syringe using clean technique.

3. POTABLE WATER RINSE:

- For all other devices, a sterile water rinse is recommended when practical. Otherwise, potable tap water rinse is acceptable.
- When using potable water for rinsing, the user should be aware of the increased risk of recontaminating the device or medical equipment with microorganisms, which may be present in potable water supplies.
- Water treatment systems, such as softeners or deionizers, may add microorganisms to the treated water to the extent that microbial content of the water at the point of use could exceed that of the pretreated drinking water. To ensure proper water quality, adherence to maintenance of the water treatment system(s) is recommended.
- The use of a bacterial retentive (0.2 micron) filter system may eliminate or greatly reduce the amount of these water borne bacteria from the potable water source. Contact the manufacturer of the filter or UV system for instructions on preventative maintenance and periodic replacement of the filter to avoid colonization or formation of biofilms in the filter.
- A device that is not completely dried provides an ideal situation for rapid colonization of bacteria. As these waterborne bacteria are highly resistant to drying, rapid drying will avoid possible colonization but may not result in a device free from these bacteria. A final rinse using a 70% isopropyl alcohol solution can be used to speed the drying process and reduce the numbers of any organism present as a result of rinsing with potable water.

C. Special Instructions for Transesophageal Echocardiography (TEE) probe reprocessing: As with all devices, carefully follow all probe manufacturer recommendations such as use of a sterile protective sheath when performing TEE. Soaking for a minimum of 5 minutes in CIDEX® OPA Solution is required for high level disinfection (HLD). Excessive soaking of the probes (e.g., longer than an hour) during HLD and/or not rinsing three times with a fresh quantity of water each time as described in Part B, may result in residual CIDEX® OPA Solution remaining on the device, the use of which may cause staining, irritation or chemical burns of the mouth, throat, oesophagus and stomach.

D. Reusage for Disinfection: CIDEX® OPA Solution has demonstrated efficacy in the presence of organic soil contamination and microbiological burden during reuse. The *ortho*-phthalaldehyde concentration of CIDEX® OPA Solution during its use-life must be verified by the CIDEX® OPA Solution Test Strips prior to each use, to determine that the MEC of 0.3% is present. CIDEX® OPA Solution may be used and reused within the limitations indicated above for up to 14 days. CIDEX® OPA Solution must be discarded after 14 days, even if the CIDEX® OPA Solution Test Strip indicates a concentration above the MEC.

MONITORING OF DISINFECTANT:

- During reuse, it is recommended that the CIDEX® OPA Solution be tested with the CIDEX® OPA Solution Test Strips prior to each use. This is to ensure that the appropriate concentration of *ortho*-phthalaldehyde is present.
- During the usage of CIDEX® OPA Solution as a high level disinfectant, it is recommended that a thermometer and timer be utilized to ensure that the optimum conditions are met.
- Visually inspect the solution during the reuse life for the presence of precipitates, which may result from the use of hard water. Discard solution if precipitation occurs.

POST-PROCESSING HANDLING AND STORAGE OF REUSABLE DEVICES:

Disinfected reusable devices are either to be immediately used, or stored in a manner to minimize recontamination. Refer to the reusable device manufacturer's labeling for additional storage and/or handling instructions.

STORAGE CONDITIONS AND EXPIRATION DATE

1. CIDEX® OPA Solution should be stored in its original sealed container at controlled room temperature 15 - 30°C (59 - 86°F) in a well-ventilated, low-traffic area.
2. Once opened, the unused portion of the solution may be stored in the original container for up to 75 days until used.
3. The expiration date of the CIDEX® OPA Solution is found on the immediate container.

EMERGENCY AND TECHNICAL PRODUCT INFORMATION

For further hazard information please refer to the Material Safety Data Sheet.

Emergency, safety, or technical information about CIDEX® OPA Solution can be obtained by contacting your local Advanced Sterilization Products sales representative.

USER TRAINING

The user should be adequately trained in the decontamination and disinfection of medical devices and the handling of liquid chemical disinfectants. Additional information about CIDEX® OPA Solution can be obtained by contacting your local Advanced Sterilization Products sales representative.

DISINFECTANT/CONTAINER DISPOSAL INFORMATION

Disinfectant Disposal: Check local disposal regulations. Glycine (free base) may be used as a neutralizer for CIDEX® OPA Solution prior to disposal. A minimum of 25 grams of glycine (free base) should be used to neutralize 3.78 litres of CIDEX® OPA Solution. The minimum recommended neutralization time is one hour. Discard residual solution into drain. Flush drain thoroughly with water.

Container Disposal: Do not reuse empty container. Rinse and dispose per hospital policy.

HOW SUPPLIED

Product Code	Description	Case Contains
20391	One Gallon (3.785L)	4 x 3.785L/case Container
20392	CIDEX® OPA Solution Test Strips	60 strips/btl; 2 btls/case
20393	CIDEX® OPA Solution Test Strips	15 strips/btl; 2 btls/case

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©ASP 2019
ASP International GmbH
EC REP **Johnson & Johnson MEDICAL GmbH**