

Simple Convenience.

TYVEK® POUCHES AND ROLLS FROM ASP. Tyvek® Pouches and Rolls with STERRAD™ Chemical Indicator are intended to be used to enclose medical devices that are to be terminally sterilized in the STERRAD™ Sterilization Systems and Cycles, in a single or double pouch configuration, and to indicate, by color change, that the pouch has been exposed to sterilant.



FEATURES AND BENEFITS

- The pouches and rolls are printed with a chemical indicator bar which is a process indicator (ISO 11140-1:2014 [Type I]) that changes from red to the color indicated on the comparator bar included on the shelf pack (or lighter) when exposed to hydrogen peroxide vapor during processing in the STERRAD™ Sterilization Systems.

Tyvek® Pouches and Rolls with STERRAD™ CI

TYVEK® SELF-SEAL POUCHES

CODE	PRODUCT	DESCRIPTION	PACKAGING
12320	Tyvek® Pouch with STERRAD™ Chemical Indicator	Self Seal, 75mm x 200mm	250/case
12326	Tyvek® Pouch with STERRAD™ Chemical Indicator	Self Seal, 100mm x 260mm	250/case
12332	Tyvek® Pouch with STERRAD™ Chemical Indicator	Self Seal, 150mm x 320mm	250/case
12335	Tyvek® Pouch with STERRAD™ Chemical Indicator	Self Seal, 100mm x 350mm	250/case
12340	Tyvek® Pouch with STERRAD™ Chemical Indicator	Self Seal, 200mm x 400mm	250/case
12342	Tyvek® Pouch with STERRAD™ Chemical Indicator	Self Seal, 150mm x 420mm	250/case
12348	Tyvek® Pouch with STERRAD™ Chemical Indicator	Self Seal, 250mm x 480mm	250/case
12356	Tyvek® Pouch with STERRAD™ Chemical Indicator	Self Seal, 320mm x 560mm	250/case



Tyvek® Self-Seal Pouches

TYVEK® HEAT-SEAL POUCHES

CODE	PRODUCT	DESCRIPTION	PACKAGING
12521	Tyvek® Pouch with STERRAD™ Chemical Indicator	Heat Seal, 75mm x 200mm	250/case
12526	Tyvek® Pouch with STERRAD™ Chemical Indicator	Heat Seal, 100mm x 260mm	250/case
12532	Tyvek® Pouch with STERRAD™ Chemical Indicator	Heat Seal, 150mm x 320mm	250/case
12541	Tyvek® Pouch with STERRAD™ Chemical Indicator	Heat Seal, 200mm x 400mm	250/case
12543	Tyvek® Pouch with STERRAD™ Chemical Indicator	Heat Seal, 200mm x 500mm	250/case
12544	Tyvek® Pouch with STERRAD™ Chemical Indicator	Heat Seal, 200mm x 650mm	250/case
12548	Tyvek® Pouch with STERRAD™ Chemical Indicator	Heat Seal, 250mm x 480mm	250/case



Tyvek® Heat-Seal Pouches

TYVEK® HEAT-SEAL ROLLS

CODE	PRODUCT	DESCRIPTION	PACKAGING
12407	Tyvek® Roll with STERRAD™ Chemical Indicator	75mm x 70m	Each
12410	Tyvek® Roll with STERRAD™ Chemical Indicator	100mm x 70m	Each
12415	Tyvek® Roll with STERRAD™ Chemical Indicator	150mm x 70m	Each
12420	Tyvek® Roll with STERRAD™ Chemical Indicator	200mm x 70m	Each
12425	Tyvek® Roll with STERRAD™ Chemical Indicator	250mm x 70m	Each
12435	Tyvek® Roll with STERRAD™ Chemical Indicator	350mm x 70m	Each
12450	Tyvek® Roll with STERRAD™ Chemical Indicator	500mm x 70m	Each



Tyvek® Heat-Seal Rolls

For more information contact your ASP representative or email: AUSNZ_ASP@asp.com

Advanced Sterilization Products Australia & ASP New Zealand
495 Blackburn Rd, Mount Waverley, Vic. 3149