



CIDEX™ OPA



EXPERIENCE THE BENEFITS

CIDEX™ OPA Solution has been cleared as high-level disinfectant for use with the most widely used endoscopes

- ➔ Multiple studies have shown the efficacy of **CIDEX™ OPA Solution** against bacteria, fungi, and viruses^{1,2}
- ➔ **CIDEX™ OPA** has been directly tested against HUMAN CORONAVIRUS, and has been demonstrated to be effective
- ➔ Low vapor pressure for minimal inhalation exposure risk
- ➔ Specially designed test strips make it easy to check the minimum effective concentration (MEC) of **CIDEX™ OPA Solution**

References

1. Akamatsu T, Minemoto M, Uyeda M. Evaluation of the antimicrobial activity and materials compatibility of orthophthalaldehyde as a high-level disinfectant. J Int Med Res. 2005;33:178-187.
2. Rutala WA, Weber DJ. New disinfection and sterilization methods. Emerg Infect Dis. 2001;7:348-353.



PROVEN. TRUSTED. SAFE.

CIDEX™ OPA

ORTHO-PHTHALALDEHYDE SOLUTION FOR HIGH-LEVEL DISINFECTION OF THE MOST WIDELY USED ENDOSCOPES

FEATURES	BENEFITS
GLUTARALDEHYDE FREE	In-use solution has low odor Effective against glutaraldehyde-resistant mycobacteria
EFFECTIVE	Effective against a wide array of microorganisms
FAST ACTING	Rapid disinfection improves productivity, allowing more endoscopes to be processed in less time
LONG LASTING UP TO 14 DAYS REUSE LIFE	Long lasting efficacy allows reprocessing of more devices per liter than with glutaraldehyde ²
EASE OF USE	CIDEX™ OPA Solution can be used straight from the bottle for extra efficiency. Requires no activation or mixing
DISPOSAL	Discard down healthcare facility drains in accordance with local regulations



Materials Compatibility

CIDEX™ OPA SOLUTION HAS BEEN TESTED WITH A WIDE RANGE OF MATERIALS COMMON TO FLEXIBLE ENDOSCOPES INCLUDING:

Aluminum, brass, copper, stainless steel, polyurethane, natural rubber latex, silicone rubber.

Refer to instructions for use for complete details.

Near neutral pH level ensures compatibility with endoscopic instruments.



ORDERING INFORMATION

REORDER NO.	DESCRIPTION	CASE CONTENTS
SCXOPA500003	CIDEX™ OPA 5 L	4 x 5L bottles
20392	CIDEX™ OPA Solution Test Strips	2 bottles (60 strips/ea.)
20393	CIDEX™ OPA Solution Test Strips	2 bottles (15 strips/ea.)

Please read and follow the CIDEX™ OPA Solution instructions for use prior to using this product for complete usage information, including contraindications and warnings.

FOR MORE INFORMATION PLEASE CONTACT YOUR LOCAL ASP REPRESENTATIVE OR VISIT ASP.COM

asp.com

ASP Advanced Sterilization Products

ASP International GmbH, Im Majorenacker 10,
Schaffhausen Switzerland
©ASP 2024. All Rights Reserved.

SCXOPA500003 | Manufactured & Packed in India
By: ALLIANCE FORMULATIONS 30B, 31-32, EPIP, Phase-I,
Jharmajri, Baddi, Solan, Himachal Pradesh, 173205

ASP International GmbH, Zweigniederlassung Zug
Gubelstrasse 34, 6300 Zug, Switzerland

EC REP ASP, The Netherlands BV
BIC 1, 5657 BX Eindhoven, The Netherlands