

YOUR SOLUTION FOR THE HIGH-LEVEL DISINFECTION AND CLEANING OF REUSABLE SEMI-CRITICAL MEDICAL DEVICES SUCH AS ENDOSCOPIC, RESPIRATORY THERAPY AND ANESTHESIA EQUIPMENT

- CIDEX™ OPA™ Solution is a high-level disinfectant for reprocessing reusable heat sensitive semi-critical medical devices
- CIDEZYME™ Enzymatic Detergent Solution is a mild enzyme-based pre-soak and cleaner designed to effectively remove organic soil
- CIDEX™ Instruments Tray Systems are used to contain cleaning and disinfecting fluid for the pre-cleaning and disinfecting of reusable medical instruments



HIGH-LEVEL DISINFECTANTS

Reorder#	Description	Packaging
SCXOPA500001 ¹	CIDEX™ OPA Solution	4 x 5L bottles
20392 ²	CIDEX™ OPA Solution Test Strips	2 bottles x 60 strips
20393 ²	CIDEX™ OPA Solution Test Strips	2 bottles x 15 strips
21396	CIDEX™ OPA Concentrate Solution	2 x 5L containers



ENZYMATIC DETERGENTS

Reorder#	Description	Packaging
2258	CIDEZYME™	12 x 1L
2260	CIDEZYME™	4 x 5 L



INSTRUMENT SOAKING TRAY SYSTEMS²

Reorder#	Description	L - W - H	Capacity
82010	Instrument Tray System w/ Perforated Inner Liner and Lid	344mm x 177mm x 132mm 13"x 7"x 5"	3.8 L
82016	Instrument Tray System w/ Perforated Inner Liner and Lid	503mm x 186mm x 132mm 20"x 7"x 5"	8 L
82027	Instrument Tray System w/ Perforated Inner Liner and Lid	740mm x 220mm x 132mm 29"x 8.5"x 5"	14 L
82032	Instrument Tray System w/ Lid and 19" (diameter) Separate Protective Compartment	449mm diameter by 222mm 19" diameter by 8.5"	17.4 L
82076	Instrument Tray System w/ Perforated Inner Liner and Lid	594mm x 381mm x 141mm 23"x 15"x 5.5"	19.7 L

Protecting patients during their most critical moments.™

asp.com

ASP International GmbH, Zweigniederlassung Zug | Gubelstrasse 34, 6300 Zug, Switzerland

²ADVANCED STERILIZATION PRODUCTS, INC. 33 Technology Drive, Irvine CA, 92618, USA

ASP Advanced Sterilization Products

ASP International GmbH, Zug Branch
Bahnhofstrasse 2, Zug 6300, Switzerland
©ASP 2022. All Rights Reserved.

EC REP ASP, The Netherlands BV
BIC 1, 5657 BX Eindhoven, The Netherlands 0123



¹Manufactured & Packed in India By: ALLIANCE FORMULATIONS 30B, 31-32, EPIP, Phase-I, Jharmajri, Baddi, Solan, Himachal Pradesh, 173205

CIDEX™ OPA and CIDEX™ OPA-C are not sterilants. Ensure adequate ventilation and wear the appropriate PPE. Please refer to the CIDEX™ OPA and CIDEX™ OPA-C instructions for use for additional information.

Call Center/ Phone: +27800064514
Email: southafrica@asp.com

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, warnings and precautions. Capitalized product names are trademarks of ASP Global Manufacturing GmbH.

AD-210014-01-CT_C-MDR