

THE ULTIMATE IN INSTRUMENT PROTECTION

Protecting your valuable, delicate surgical instruments and devices during sterilization, storage, and transport is vital. That's why Advanced Sterilization Products offers a full line of trays, silicone mats and instrument holders to secure and protect your instruments sets. APTIMAX™ Instrument Trays feature a rugged, patented material that is compatible with most commonly used sterilization systems. It's easy to see why APTIMAX™ Instrument Trays and Accessories offer the ultimate in protection for your valuable instruments.



FEATURES AND BENEFITS

Protection

- Durable material stands up to washer and sterilization processes

Standardization

- Fully validated for use in all STERRAD™ Systems, pre-vacuum steam and ethylene oxide sterilization
- Full line of trays, mats and instrument holders allow you to customize instrument sets

Convenience

- Lightweight design and snug-fitting lid for ease in handling
- Stackable design to optimize space efficiency for storage

APTIMAX™ INSTRUMENT TRAYS & ACCESSORIES

APTIMAX™ INSTRUMENT TRAYS


Code	Product description	Dimensions	Packaging
13826	APTIMAX™ Instrument Tray	300 mm x 100 mm x 50 mm	2/case
13827	APTIMAX™ Instrument Tray	450 mm x 100 mm x 50 mm	2/case
13828	APTIMAX™ Instrument Tray	600 mm x 100 mm x 50 mm	2/case
13829	APTIMAX™ Instrument Tray	300 mm x 200 mm x 50 mm	2/case
13830	APTIMAX™ Instrument Tray	450 mm x 200 mm x 50 mm	2/case
13831	APTIMAX™ Instrument Tray	600 mm x 200 mm x 50 mm	2/case
13832	APTIMAX™ Instrument Tray	180 mm x 75 mm x 30 mm	2/case
13833	APTIMAX™ Instrument Tray	250 mm x 150 mm x 30 mm	2/case
13834	APTIMAX™ Instrument Tray	380 mm x 150 mm x 30 mm	2/case
13835	APTIMAX™ Instrument Tray	380 mm x 250 mm x 30 mm	2/case
13836	APTIMAX™ Instrument Tray	278 mm x 270 mm x 100 mm	2/case
13837	APTIMAX™ Instrument Tray	576 mm x 270 mm x 100 mm	2/case
13838	APTIMAX™ Instrument Tray	596 mm x 321 mm x 100 mm	1/case



APTIMAX™ Instrument Tray

STERRAD™ INSTRUMENT TRAY MATS


Code	Product description	Dimensions	Packaging
99205	APTIMAX™ Instrument Tray Mat	165 mm x 64 mm fits tray 13832	2/case
99206	APTIMAX™ Instrument Tray Mat	229 mm x 140 mm fits tray 13833	2/case
99207	APTIMAX™ Instrument Tray Mat	356 mm x 140 mm fits tray 13834	2/case
99208	APTIMAX™ Instrument Tray Mat	356 mm x 241 mm fits tray 13835	2/case
99209	APTIMAX™ Instrument Tray Mat	279 mm x 191 mm fits tray 13829	2/case
99210	APTIMAX™ Instrument Tray Mat	432 mm x 191 mm fits tray 13830	2/case
99211	APTIMAX™ Instrument Tray Mat	584 mm x 191 mm fits tray 13831	2/case
99212	APTIMAX™ Instrument Tray Mat	254 mm x 254 mm fits tray 13836	2/case
99213	APTIMAX™ Instrument Tray Mat	546 mm x 254 mm fits tray 13837	2/case



APTIMAX™ Instrument Trays & Mats

APTIMAX™ INSTRUMENT TRAY HOLDERS

Code	Product description	Dimensions	Packaging
99400	APTIMAX™ Instrument Tray Holder	10 mm x 19 mm x 10 mm	4/case
99401	APTIMAX™ Instrument Tray Holder	19 mm x 10 mm x 19 mm	4/case
99402	APTIMAX™ Instrument Tray Holder	5 mm x 19 mm x 5 mm	4/case
99403	APTIMAX™ Instrument Tray Holder	19 mm x 5 mm x 19 mm	4/case
99404	APTIMAX™ Instrument Tray Holder	10 mm x 10 mm x 10 mm	4/case
99405	APTIMAX™ Instrument Tray Holder	15 mm x 15 mm x 15 mm	4/case
99406	APTIMAX™ Instrument Tray Holder	5 mm x 5 mm x 5 mm	4/case
99407	APTIMAX™ Holding Bar		4/case



APTIMAX™ Instrument Tray Holder

Protecting patients during their most critical moments.™