

Innovation to advance the Quality of Patient Care

Discover the latest suite of technologies from ASP



Enhance Compliance Automatically

Discover an unprecedented level of interconnectivity that facilitates productivity and enhances compliance.*

ASP technologies now work together to help prevent mistakes before they happen. On-screen, step-by-step instructions with automatic checkpoints help enable compliance with Instructions for Use (IFU) and hospital policy.

- Biological indicator (BI) usage reminders
- Cycle selection and load placement instructions
- Test BIs linked with corresponding control BI
- Automatic and immediate reconciliation of sterilizer cycle and BI results
- Automated record keeping
- Records accessible and audit-ready



ASP ACCESS™ Technology

Smart information-sharing technology that enhances compliance and easily integrates with instrument tracking systems (ITS) and hospital networks



STERRAD VELOCITY™ Biological Indicator System

Results in 15 minutes enable healthcare facilities to reduce the risk of releasing instruments prior to BI confirmation



STERRAD™ Systems with ALLClear™ Technology

ALLClear™ Technology minimizes workflow interruptions and user errors to drive compliance

*ASP ACCESS™ Technology, STERRAD™ Systems with ALLClear™ Technology, and the STERRAD VELOCITY™ Biological Indicator System have features that may enhance compliance, including BI reminders per hospital policy, and on-screen reinforcement of user training.

STERRAD™ Systems with ALLClear™ Technology JUST LOAD AND GO

Trusted sterilization from the only system with automated reconciliation of sterilization cycle and BI records



Now featuring ALLClear™ Technology, the STERRAD™ 100NX System and the STERRAD NX™ System are products designed to maximize efficiency and compliance* in the CSSD.

01 Reduce workflow interruptions

- Fast, accurate detection and correction of load and system issues



System Check

System diagnostics process completes prior to the sterilization cycle



Moisture Detection & Correction

Automatically detects and removes moisture



Instrument Warming

Plasma transfers energy into the load, removing gas and warming instruments

02 Enhance compliance automatically*

- Quality control that enhances compliance with hospital policy and IFUs
- Minimize the need for manual record keeping

03 Designed with simplicity in mind, to reduce the potential for human error

- Enhanced user interface enables fast and easy staff training



*STERRAD™ Systems with ALLClear™ Technology have features that may enhance compliance, such as BI reminders per hospital policy, and on-screen reinforcement of user training.

ASPTM





STERRAD VELOCITY™

Biological Indicator / Process Challenge Device

RESULTS IN 15 MINUTES*

Smart BI/PCD technology that is fully integrated with all the latest features of STERRAD™ System sterilizers

- The first all-in-one Process Challenge Device designed by ASP specifically for all STERRAD™ System Cycles to minimize false positives
- Enables healthcare facilities to reduce the risk of releasing instruments prior to BI confirmation
- Provides the protection of the only rapid read Process Challenge Device (PCD) for STERRAD™ Systems that meets the elevated standard of AAMI recommended guidelines (ST58 9.5.4.3)



*15 or 30 minutes to result dependent on the software version on the STERRAD VELOCITY™ Reader



Engineered to be easy to use


- Intuitive on-screen guidance enable correct operation, reducing the incidence of user error



- Reduces the need for manual record keeping through automated correlation of BI results and sterilization cycle records¹

- Communicates across STERRAD™ Systems, ITS, and hospital networks through connection with ASP ACCESS™ Technology

Enhanced compliance* through connectivity



*The STERRAD VELOCITY™ Biological Indicator System has features that may enhance compliance, including BI reminders per hospital policy, as well as improved record handling and record keeping. ¹Automated integration of information requires connection to ASP ACCESS™ Technology.

ASP ACCESS™ Technology

CONNECT. CORRELATE. COMMUNICATE.

ASP ACCESS™ Technology facilitates communication between STERRAD™ Systems, STERRAD VELOCITY™, ITS, and hospital networks

Delivers smart information-sharing technology designed to enhance compliance* and provide unique workflow insights

01

Automatically and immediately communicates important sterilizer cycle results and system status changes to connected devices

→ Cloud technology enables staff to receive alerts via email and text

02

Automatically and immediately communicates important sterilizer cycle results and system status changes to connected devices

→ Seamless collection of vital reprocessing information reduces the risk of human error by simplifying the documentation process

→ Cloud technology enables immediate access to sterilization records for multiple departments and hospital sites

ITS

SISTEMA DI
TRACCIABILITÀ
DEGLI STRUMENTI



*ASP ACCESS™ Technology has features that may enhance compliance, including automatic reconciliation of sterilization records and communication with hospital networks and instrument tracking systems, to improve record handling and record keeping.

Protecting Patients in Their Most Critical Moments



asp.com

ASP Advanced Sterilization Products

Sede Operativa
Advanced Sterilization Products Italia Srl
Via Vaccareccia 39 A - 00071 Pomezia (Roma)
©ASP 2021. Tutti i diritti riservati.



ADVANCED STERILIZATION PRODUCTS, INC.
33 Technology Drive, Irvine CA 92618, USA

EC REP

ASP, The Netherlands BV
BIC 1, 5657 BX Eindhoven, The Netherlands

CE
0123

Please read and follow the Instructions for Use (IFU). Important information: Prior to use, refer to the complete instructions for use supplied with the device(s) for proper use, indications, contraindications, warnings and precautions.
Capitalized product names and ALLClear™ are trademarks of ASP Global Manufacturing, GmbH.
For more information, please contact your local ASP representative or visit asp.com

AD-170035-01-CT_D-MDR